


CXA000 12' Swing Sampler w/Pole and Bottle

The CXA000 Swing Sampler is designed for collecting liquid samples from a horizontal flowing stream such as a sewer. It allows collection from various angles, including 90°. Different sizes of bottles/containers may be attached to the bottle holder, allowing for a variety of applications. A 960 ml bottle is provided with the Swing Sampler and Global Water also offers a 500 ml bottle and corresponding snapper ring. The pole can be extended up to 12 feet (3.66 m).



Assembly Instructions

To connect the pole to the bottle holder:

1. Locate part B (on the pole) and part C (on the bottle holder).
2. Turn part B clockwise into part C until it is snug.
3. Use the supplied Allen wrench to tighten the Allen screw located in the inner pole screw threads of the bottle holder.

Assembly is now complete